

TL-232P

Trial lens set



MAGNON

TL-232P Trial lens set

The MAGNON TL-232P is double curve trial lens set which offers a standard selection of lens powers to meet practitioner's need.

FEATURES

- Double curve lenses of 22mm in diameter
- Color coded plastic rims of 38mm in diameter
- Removable plastic inner tray with space for trial frame
- Aluminum carrying case with combination locks

SPECIFICATIONS

Sphere Cylinder Prism Accessories

- 34 pairs each of plus and minus
19 pairs each of plus and minus
10 pieces
10 pieces
- : Occluder 2
 - : Pin hole (1mm) 1
 - : Pin hole (1.5mm) 1
 - ⊖ : Silt (1.5mm) 1
 - R : Red lens 1
 - G : Green lens 1
 - F : Frosted lens 1
 - N : 0.0D lens 1
 - M : White Maddox 1

- Spheres	- Cylinders	+ Cylinders	+ Spheres
20	6.00	6.00	20
18	5.00	5.00	18
16	4.50	4.50	16
15	4.00	4.00	15
14	3.50	3.50	14
13	3.25	3.25	13
12	3.00	3.00	12
11	2.75	2.75	11
10	2.50	2.50	10
9.00	2.25	2.25	9.00
8.00	2.00	2.00	8.00
7.00	1.75	1.75	7.00
6.50	1.50	1.50	6.50
6.00	1.25	1.25	6.00
5.50	1.00	1.00	5.50
5.00	0.75	0.75	5.00
4.50	0.50	0.50	4.50
4.00	0.25	0.25	4.00
3.75	0.12	0.12	3.75
3.50	Accessories	Prisms	3.50
3.25	●	1°	3.25
3.00	●	2°	3.00
2.75	○	3°	2.75
2.50	○	4°	2.50
2.25	⊖	5°	2.25
2.00	R	6°	2.00
1.75	G	7°	1.75
1.50	F	8°	1.50
1.25	N	9°	1.25
1.00	M	10°	1.00
0.75			0.75
0.50			0.50
0.25			0.25
0.12			0.12

*Specifications and design are subject to change for improvement without notice.

H. OGINO & CO., LTD.

Yokohama Port P. O.Box 170, Yokohama 231-8691, Japan
Tel: (045)662-5261 Fax: (045)651-5467
E-mail: magnon@h-ogino.co.jp URL: www.h-ogino.co.jp



OMMALITE
OMMALITE LENSES & OPTICAL INSTRUMENTS

Μεταμορφώσεως & Χίου 6-8, 183 45, Μοσχάτο
τηλ: 210 48 19 098 / fax: 210 48 19 946
Γαμβέτα 195B, 542 48, Θεσσαλονίκη
τηλ: 2310 30 04 67 / fax: 2310 31 90 52

www.ommalite.gr
info@ommalite.gr

01/06/13A Printed in Japan.