





www.yuratek.com

OPHTHALMIC PRODUCTS

For more than 14 years Yuratek has been developing products according to its experiences in ophthalmology and optometry industries all around the world and exporting products to many countries in Middle East, Europe, Asia and Africa. Our primary goal is to our customers happy with the quality and customizable options...

BAHA II



BAHAIII



YRT 2000



YRT 1000



TABLES



BAHA II



10-250-0>

BAHA II offers a comfortable configuration as well as flexible expandability.

BAHA II includes a swivel-sliding instrument table for two instruments and process a drawer element offering room for trial lens set and other accessories.

The hinged drawer element offers the best working position during refraction every time.

The easily accessible control panel serves for simple and fast operation of BAHA series.





The double column that is added to the unit table for the motorized up and down movement of the phoropter arm adds speed and comfort to the doctor examination.



- * The wheel chair version is a factory option.
- Unit stabilite is guaranteed by extra weight to the base plate.
 Designed according to ISO standards for wheel chair configuration.



Left or right sided table installation feature allows you to easily adjust the unit's table during the installation.



Foldable Footrest

Standard Features

- · Left and right sided selection feature during installation
- Magnetic lock mechanism (can be locked in every position)
- · Sliding and rotating table top for two instruments
- Chinrest preparation compatible with almost all brand names
- · Chair up-down movement
- · Power outlets including AUX and hidden cable channels
- · Patient knee protection switch
- Easy to control keypad (changeable left or right side)
- Chart projector column
- · Chart projector holder and plate
- · Dimmable halogen head light
- · Ergonomic chair with footrest, armrest
- · Wide range of color and tissue
- All units in the BAHA series are equipped with an electrically operated patient's chair.
- Push-button for descending chair. A fast double pulsation, pressing twice make descend chair automatically to the response position.



Controls available for:

Power On-Off

Reading Light (Led Lamp)

Dimming of Slit Lamp

Dimming of Head Lamp

Chair Up-Down, Automatic Landing "0" Position

Table Up-Down

Aux On-Off

Projector On-Off

Motorized Phoropter Arm



The projector is mounted on the aluminium column.

BAHA II





Tabletop Up - Down Movement



Triple Drawer



Horizontal Chair Movement Manual or Motorized



Rotating Chair 360°



Inclinable Chair max your patient feels like home 45° or full

12-400

Motorized Sliding Phoropter Arm



12-300

Manual Sliding Phoropter Arm



12-200

Motorized Phoropter Arm



12-100

Manual Phoropter Arm



12-500

Manual Phoropter Arm with Wall Adapter



10-500

LCD Monitor Support



00-200


Foot Pedal for Electrical Chair





00-300


Electro-Brakes are Controlled for a Foot Switch











BAHA III



For highly equipped examination rooms BAHA III series offer a complete solution with rotating almost 355 degrees tables and optional table up and down movement.

Standard Features


- Tabletop for 3 instruments
- Both right and left units are available with table top elevation
- Combination of two smooth and fast acting movements of table top (both rotation and sliding 355°)
- Which aporete with electro-magnetic brake to place each instrument in the correct position (every angle)
- · Dimming of slit lamp control panel
- · Chair for without armrest
- · Standard trail lens drawer
- Wide range of color and tissue selection for both unit body and chair cover
- Head lamp is upper plate aluminum port of the halogen lamp of 230V/50W and dimmable control main panel
- Extruded aluminum with side slots to attach the projector and phoropter arms, both height adjustable head lamp projector suport
- · Control panel with an easy access
- Powered elevation of tabletop (0-6/12V and 230 V Ac)




Left or Right sided table installation feature allows you to easily adjust the unit's table during the installation.

Optional Accessories

- Right and left wheel chair access version available (factory option)
- · Up-down movement for instrument table top
- Wheelchair access possibility (factory option)
- Single trial drawer or a three drawer
- Inclinable chair backrest 45° to full
- Electro-brakes whics are controlled by a foot switch
- Horizontal chair movement manual or motorized
- Rotating chair 360°
- LCD monitor support
- Ophtalmoscope retinascope holder
- Choice out of phoropter arm with regulation for different weights (4-7 kg.) and adjustable in force Horizontal & vertical position and direction. (code 12-100)





YRT 2000



Ergonomics and economics are keywords in the YRT-2000 opthalmic unit

Standard Features

- · Easy to operate control panel
- Dimming of head light
- Dimming of slit lamp
- Projector on/off
- One auxiliory instrument
- · Motorized chair with arm rest and footrest
- · Mechanical lock for sliding table
- · One trail lens drawer
- · Projector colum and support table
- Unit in the YRT 2000 is equipped with an electrically operated patient's chair
- Push-button for descending chair. A fast double pulsation, pressing twice make descend chair automatically to the respose position
- Mechanical lock for sliding table









Optional Accessories

- Double drawer
- Electromagnetic lock
- · Choice out of four type phoropter arm
- Electromagnetic brake on/off
- Chair rotation 360°
- Inclaniable chair 45° or full
- Left side assembly possibility (factory option)
- LCD monitor support
- Oftolmascope retinascope holder








TABLES




| Specifications | MT-30 | |
|------------------|-----------------------------------|--|
| Input Voltage | 240 VAc (Optional 100-240 VAc) | |
| Power Aux Outlet | 2 | |
| Maximum Load | 100 Kg | |
| Stroke | 250 mm | |





| Specifications | MT-20 |
|------------------|-----------------------------------|
| Input Voltage | 240 VAc (Optional 100-240 VAc) |
| Power Aux Outlet | 2 |
| Maximum Load | 75 Kg |
| Stroke | 250 mm |









| Specifications | MT-10/MT-10-S |
|------------------|-----------------------------------|
| Input Voltage | 240 VAc (Optional 100-240 VAc) |
| Power Aux Outlet | 1 |
| Maximum Load | 75 Kg |
| Stroke | 250 mm |









CHAIR











Fixed backrest

Fixed backrest


Tiltable backrest 45° or full

SPECIFICATIONS 11-500 11-700 11-750 Fixed Backrest Yes Yes No No Tiltable Backrest No 45° or full Foldable Armrest Yes Yes Yes Foldable Footrest Yes Yes Yes Optional Optional Optional Horizontal Seat Movement

Motorized Chair







Fixed backrest
Electrical foot pedal control
Up-Down movement 180 mm
Max, load capacity 160 kg
Power input 230 V Ac/50 Hz
Optional 100-230 V Ac/50 Hz



Tiltable backrest 45° or full Electrical foot pedal control Up-Down movement 180 mm Max. load capacity 160 kg Power input 230 V Ac/50 Hz Optional 100-230 V Ac/50 Hz


Optional Headrest



| SPECIFICATIONS | BAHA II | BAHA III | YRT 2000 | YRT 1000 |
|--|------------------------------------|------------------------------------|--|--------------|
| Tabletop | 2 Instruments | 3 Instruments | 2 Instruments | 1 Instrument |
| Sliding tabletop | Rotatable sliding tabletop | <u> </u> | Rotatable sliding tabletop | |
| Maximun load tabletop | 40 kg (20 + 20 kg) | Each table top : 20 kg | 40 kg (20 + 20 kg) | 20 kg |
| Max. rotation angle swivel tabletop | 90° | 355° | 90° | 90° |
| Wheelchair access | Available | Available | | Usable |
| Seat rotation | 0 to 360° | 0 to 360° | 0 to 360° | 7.4 |
| Maximun load for chair | 175 kg | 175 kg | 175 kg | |
| Electro-brake for locking the tabletop | Can be locked in every position | Can be locked in every position | Can be locked in every position (optional) | |
| Power head lamp | 220 V Ac/50W Halogen (Dimmable) | 220 V Ac/50W Halogen (Dimmable) | 220 V Ac/50W Halogen (Dimmable) | Tes |
| Stroke elevated tabletop | 780-950 mm | 800-980 mm | fixed-880 mm | 770-950 mm |
| Trail lens drawer | 1 or 3 | 1 or 3 | 1 or 2 | |
| Inside dimensions trail lens drawer | 530(w) x 360 (d) x 70(h) | 530(w) x 360 (d) x 70(h) | 530(w) x 350 (d) x 70(h) | |
| Total power supplay | 220-240V / 50-60 Hz | | | |
| Power near vision lamp | LED (Optional) | | | |
| Safety stop | YES | | | |
| Power on tabletop | 0-6V-12V-16 V/5A 220VAc-50/60 Hz | | | |
| Maximum power consumption | 900 W | | | |
| Aux voltage | Max 240 VAc / 5A | | | |
| Electrical classification | Class I, Type B | | | |
| Electric protection degree | | IP | 21 | |

Wood Box Dimensions for BAHA II - BAHA III - YRT 2000

| BAHA II | 250 kg. |
|----------|---------|
| BAHA III | 275 kg. |
| YRT 2000 | 225 kg. |
| YRT 1000 | 125 kg. |





COMPLY TO THE DIRECTIVES

93/ 42/ EEC Medical Device Directive (MMD)
2007/47/EC Revision to Medical Device Directives
2006/95/EC Low Voltage Directive
2011/65/EU EU RoHS 2 Directive

UNIT COLOURS OPTIONS



CHAIR COLOURS OPTIONS

| 001 | 013 | 016 | 003 |
|----------------|----------------|----------------|----------------|
| | | | |
| 030 | 036 | 210 | 215 |
| | | | |
| 220 | 227 | 600 | 345 |
| | 415 | 425 | |
| 410 | 410 | 425 | 590 |
| 520 | 530 | 580 | 582 |
| | | | |
| 515 | 320 | 620 | 901 |
| Sample Leather | Sample Leather | Sample Leather | Sample Leather |



lvedik Org. San. Bölgesi 1475 Sk. (eski 692. Sk.) No: 4 Ostim - ANKARA - TURKEY Tel: +90.312 395 29 96 • Fax: +90.312 395 29 97 • e-mail: info@yuratek.com



Basım Tarihi: 16 ŞUBAT 2015 Basım Yeri:

yalçınmatoaacılık wambalajsanal Erman BÜYÜKAKGÜL www.yalcinmatbaacilk.comtr T. 0312 435 0 435 - 395 21 16